


Highly Sensitive and Comprehensive Detection for Chiral Separation of Pesticide using HPLC-CD-MS

Introduction

About 30% of pesticides currently used worldwide are chiral isomers with one or more chiral centers. It is well known that the enantiomers of chiral pesticides have different biological activities, toxicity on non-targeted organisms, rate of metabolism, and biodegradation in the environment. Although it is desirable to use the more effective and safer enantiomers, most chiral pesticides are marketed and used as racemates due to the complicated process and high cost of manufacturing a single enantiomer molecule.

In general, chiral separation by HPLC, GC, and SFC has been one of the most important ways for analyzing the synthetic approach to individual enantiomers. However, HPLC still dominates chromatographic chiral analysis because of its ease of use.

In this study, we developed a chiral separation method for pesticides using HPLC with circular dichroism detector (CD) and mass spectrometer (MS). The CD detector simultaneously provides circular dichroism and UV absorbance data, and enables selective identification of the CD polarity of each enantiomer. MS spectra are used to distinguish between chiral isomer groups with the same m/z from other components including impurities. Highly sensitive and comprehensive detection with UV, CD, and MS contributes to the widely applicable and highly accurate analysis for chiral separation in the efficient development of a single enantiomer pesticide. The system and chiral separation method described here is applicable to many industries that deal with chiral compounds. In this presentation, we will present the application of HPLC to bromuconazole, a type of triazole fungicide (Figure 1).



MW 377.06

Figure 1. Structure of Bromuconazole

Experimental

Apparatus

Figure 2 details the JASCO HPLC-CD-MS system used in this experiment, together with Figure 3 that shows the schematic diagram of the system. Using a splitter with a PEEK-coated capillary tube, the flow from a chiral column is divided into two channels, one each to the CD and MS detectors. The CD detector provides CD and UV signals simultaneously. This system enables automated method scouting analysis to be performed on mobile phases and columns.




Mobile Phase Pump	PU-4180 (with DG ¹ and LPG ² unit)
Make-up Pump	PU-4185 (with DG ¹ unit)
Autosampler	AS-4050
Column Oven	CO-4065 (with 2 switching valve unit ³)
Detector 1	CD-4095
Detector 2	expression CMS (single quadrupole)
Chromatography Data System	ChromNAV Ver.2

¹ Degassing unit
² Low pressure gradient unit
³ 6 position-7port valve unit

Figure 2. The JASCO LC-CD-MS system

EXTREMA



- 1 Mobile phase
- 2 Degassing unit
- 3 Low pressure gradient unit
- 4 Mobile phase pump
- 5 Autosampler
- 6 Column oven (with 2 switching valve units)
- 7 Chiral columns
- 8 Flow splitter
- 9 Make up solvent for MS
- 10 Make up pump
- 11 Circular dichroism detector
- 12 Mass spectrometer

Figure 3. Schematic diagram of the LC-CD-MS.

Analysis Conditions for Method Scouting

Chiral Columns	1; CHIRALPAK AD-H	2; CHIRALPAK AS-H	3; CHIRALCEL OD-H
Column Size	4; CHIRALCEL OJ-H	5; CHIRALPAK IA	6; CHIRALPAK IB
Mobile Phase	4.6 mmI.D. x 150 mmL, 5 μm		
Mobile Phase Composition	A; n-Hexane (Hex)	B; Ethanol (EtOH)	C; 2-propanol (IPA)
Make-up Solvent	Ethanol / Water (95 / 5, with 0.1 % formic acid)		
Flow Rate	Mobile phase ; 1.0 mL/min, Make up solvent ; 0.2 mL/min		
Column Temperature	30° C		
CD Detection	230 nm		
MS Detection	Ionization mode ; ESI-positive		
	Measurement mode ; Scan		
	Mass range ; 100-600 m/z		
	Capillary temperature ; 250° C		
	Capillary voltage ; 160 V		
	Source gas temperature ; 350° C		
	ESI voltage ; 3500 V		
Sample	500 μg/mL racemic bromuconazole in hexane/ethanol (90/10)		
Injection Volume	10 μL		
Split Ratio	CD: MS = ca. 15 : 1		

Results and Discussion

Results of Method Scouting on Mobile Phases and Columns

Racemic bromuconazole with 2 chiral centers can be separated into 4 isomer peaks using a chiral column. Figure 4 shows the selected results of method scouting analysis for mobile phases and columns. The best separation was observed with the combination of n-hexane/ethanol (90/10) and CHIRALPAK AS-H (highlighted in yellow). This separation condition was used in the following experiments.




Figure 4. Chromatograms of bromuconazole standard obtained by method scouting analysis, (A) CD chromatograms, (B) UV chromatograms. Rows show mobile phases and columns show chiral columns.

Comprehensive Detection by CD and MS Detectors

Figure 5 shows the results of bromuconazole standard spiked with two additional pesticides, pirimicarb and triadimefon, as examples of contaminants detected by CD and MS detectors. Figure 5A shows the total ion chromatogram (TIC) from 100 to 600 m/z, extracted ion chromatogram (XIC) of m/z 378.1, CD and UV chromatograms. Figure 5B shows MS spectra of each peak described in Figure 5A. As shown in this figure, we can clearly discriminate the bromuconazole isomers and other contaminants. Bromuconazole isomers (peaks 3 ~ 6) have similarities in their spectra, and a protonated molecular ion ($[M+H]^+$) was observed at m/z 378.1 as a base peak in the mass spectrum of each isomer peak. Protonated molecular ions of pirimicarb and triadimefon are also observed at m/z 239.2 (peak 1) and m/z 294.2 (peak 2). Figure 5C shows the stopped-flow scanned spectra of CD and UV for bromuconazole standard peaks (peaks 3 ~ 6). According to the CD information in Figure 5A and Figure 5C, peaks 3 and 5 of the bromuconazole isomers are estimated to be an enantiomeric pair from their CD polarity and wavelength maxima. Peaks 4 and 6 are also the same. Furthermore, peak 2 in the CD chromatogram in Figure 5A indicates the possibility that can be each enantiomer of triadimefon can be separated.

In this way, comprehensive detection by CD and MS detectors is useful to identify the target chiral isomers group and other contaminants with chiral separation.




Figure 5. The results of a bromuconazole standard spiked with contaminants using CD and MS detectors, (A) MS, CD, and UV chromatograms, (B) MS spectra of each peak, (C) stopped-flow scanned spectra of CD and UV for bromuconazole peaks (peak 3 ~ 6). Racemic bromuconazole: 500 $\mu\text{g}/\text{mL}$, pirimicarb: 50 $\mu\text{g}/\text{mL}$, racemic triadimefon: 100 $\mu\text{g}/\text{mL}$.

Linearity, Reproducibility and Sensitivity

Figure 6 shows the overlaid chromatograms of SIM (Single Ion Monitoring) at m/z 378.1, CD, and UV of bromuconazole standards (2.5, 12.5, and 25 $\mu\text{g}/\text{mL}$ enantiomer concentration). Table 1 shows the linearity, reproducibility of retention time and peak area, and detection limit of the bromuconazole enantiomers (the concentrations are for enantiomers). As shown by these results, a wide linear dynamic range and highly sensitive detection were obtained, especially for the SIM chromatogram.




Figure 6. SIM, CD, and UV chromatograms of bromuconazole standards (2.5, 12.5, and 25 $\mu\text{g/mL}$ of enantiomer concentration).

Table 1. Linearity, reproducibility of retention time and peak area, and detection limit of bromuconazole enantiomers

	Correlation Coefficient	%RSD (n=10)		Detection Limit (S/N=3)
		tR	Peak Area	
MS-SIM (<i>m/z</i> 378.1)	0.9994 ~ 0.9997 (0.025 ~ 25 $\mu\text{g/mL}$)	0.06 ~ 0.11 ¹	1.51 ~ 1.99 ¹	0.0034 ~ 0.0057 $\mu\text{g/mL}$ ³
CD	0.9987 ~ 0.9994 (12.5 ~ 125 $\mu\text{g/mL}$)	0.06 ~ 0.14 ²	0.86 ~ 2.64 ²	0.55 ~ 5.06 $\mu\text{g/mL}$ ⁴
UV	0.9992 ~ 0.9997 (2.5 ~ 125 $\mu\text{g/mL}$)	0.06 ~ 0.11 ²	1.11 ~ 1.31 ²	0.089 ~ 0.16 $\mu\text{g/mL}$ ⁴

¹ Calculated from the results of 1.25 $\mu\text{g/mL}$.

² Calculated from the results of 62.5 $\mu\text{g/mL}$.

³ Calculated from the results of 0.125 $\mu\text{g/mL}$.

⁴ Calculated from the results of 25 $\mu\text{g/mL}$.

Conclusion

We developed a chiral separation method for pesticides by HPLC with CD and MS detection.

Detection by TIC, CD, and UV enable identification of the target chiral isomers group and other contaminants with chiral separation.

Wide dynamic range and high sensitivity for each bromuconazole enantiomer were obtained, especially in the SIM chromatogram.

This system is applicable to many industries that deal with chiral compounds.

References

B. S. Sekhon, J. Pestic. Sci, 2009, 34 (1), 1-12.

V. Pérez-Fernández, M. Á. Garcíá, M. L. Marina, J. Chromatogr. A, 2011, 1218, 6561-6582.