


### LC/MS Analysis of Aminoglycoside Antibiotics (VC-50 2D)

Aminoglycoside antibiotics have four amino groups making them highly hydrophobicity. When they are analyzed by reversed phase mode, it is necessary to add an ion pair reagent into the eluent. These antibiotics do not have UV absorption, so LC/MS method is suitable for high sensitivity analysis of them. In this application, various aminoglycoside antibiotics were analyzed using HILICpak VC-50 2D, a polymer-based HILIC mode column, in cation exchange mode with LC/MS detection. Since alkaline conditions suppress the dissociation of amino groups in the molecules, the cation exchange between the antibiotics and gel materials weakens, allowing the antibiotics to elute.



Sample : Aminoglycoside antibiotics 0.1 $\mu$ g/mL (in H<sub>2</sub>O), 20  $\mu$ L

Kanamycin sulfate  
Tobramycin sulfate  
Gentamicin sulfate  
Neomycin sulfate  
Amikacin dihydrate

Column : Shodex HILICpak VC-50 2D (2.0mmI.D. x 150mm)  
Eluent : (A); 1.5% NH<sub>3</sub> aq./ (B); CH<sub>3</sub>CN  
High pressure linear gradient ;  
(B%) 30% to 10% (0 to 5min), 10% (5 to 15min)  
Flow rate : 0.3mL/min  
Detector : ESI-MS SIM  
Column temp. : 40°C