

# Micro Weighing Funnels

Precision Meets Convenience – Effortless Sample Transfer Every Time.

## No Spills. No Loss.

The micro measuring funnel is a convenient tool for weighing samples using a microbalance. After weighing, the sample can be easily transferred into a vial by inserting the funnel's narrow stem into the vial opening. This allows solvent to be poured in smoothly, making sample transfer much easier and preventing any loss of material.



## How to Use the Micro Weighing Funnel

### Step 1.

Weigh the sample on a microbalance. Place the funnel on the balance and carefully add the sample directly into the funnel.




### Step 2.

Insert the funnel into a test tube or vial and pour solvent from one side. Then gently pour the solvent along the inner wall of the funnel to dissolve and transfer the sample without spillage.



## Ordering Information

Descriptions	Size A (mm)	Size B (mm)	Size C (mm)	Size D (mm)	Weight (g)	Cat. No.
Weighing Funnel Set (1mL, 2 mL, 10 mL)	–	–	–	–	–	1050-11000
Weighing Funnel 1 mL	28	35	5	10	Approx. 2.0	1050-11001
Weighing Funnel 2 mL	32	42	6	13.7	Approx. 4.2	1050-11002
Weighing Funnel 5 mL	35	47	8	18.5	Approx. 6.6	1050-11005
Weighing Funnel 10 mL	45	58	8	23	Approx. 10	1050-11010
Weighing Funnel 20 mL	55	70	8	29	Approx. 19	1050-11020



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement. Please note that, in the interests of continuous improvement, models or specifications are subject to change without notice. Please also note that the company name and product name appearing in this catalogue are the trademark or registered trademark of each corresponding company. In the descriptions in this catalogue, TM and R marks are not used.

Contact us or your local GL Sciences representative.

<https://www.glsciences.com/contactus/index.php>

Authorized distributor:

<https://www.glsciences.com/company/distributor.html>

For analytical use only. Information subject to change without prior notice.

 **GL Sciences** <https://www.glsciences.com>  
E-Mail [world@glsc.co.jp](mailto:world@glsc.co.jp)