

## Lorcaserin

### CHROMATOGRAPHIC CONDITIONS

Column: CHIRALPAK® IF [4.6 x 250 mm, 5 µm]  
Part No: 86325

Mobile phase: Methanol : EDA (100/0.1, v/v)


Flowrate: 1.0 mL/min

Column Oven Temperature: 25 °C

UV Detection: 230 nm

Diluent: Mobile phase

Sample Concentration: 2.0 mg/mL



### SEPARATION PERFORMANCE

$t_{R1}$ : 5.5	$N_1$ : 8687	$T_1$ : 1.5	
$t_{R2}$ : 6.5	$N_2$ : 5776	$T_2$ : 1.9	$R_s$ : 3.3

