

Rebamipide

CHROMATOGRAPHIC CONDITIONS

Column: CHIRALPAK® IE-3 [4.6 x 250 mm, 3 µm]

Part No: 85525

Mobile phase: ACN : Methanol : TFA (85/15/0.3, v/v/v)


Flowrate: 0.7 mL/min

Column Oven Temperature: 40 °C

UV Detection: 230 nm

Diluent: Methanol : ACN (1:1)

Sample Concentration: 0.2 mg/mL



SEPARATION PERFORMANCE

t_{R1} : 6.9	N_1 : 14300	T_1 : 1.3	
t_{R2} : 7.9	N_2 : 10807	T_2 : 1.3	R_s : 3.7

